

OsteoCentric®

Integrated Locking Proximal Humerus System

UnifiMI® technology resists multi-planar and multi-directional loads through Mechanical Integration.

UnifiMI® is designed to engage twice as much bone when loaded in any direction as compared to standard screws.

UnifiMI® 4.0mm screw fastener improves bone capture and engagement in good and compromised bone.



OsteoCentric®
TECHNOLOGIES

This document is intended exclusively for physicians and is not intended for laypersons. Information on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendations. Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and/or advice.

All content herein is protected by copyright, trademarks and other intellectual property rights owned by or licensed to OsteoCentric Technologies, Inc. or one of its affiliates and must not be redistributed, duplicated, or disclosed, in whole or in part, without the express written consent of OsteoCentric Technologies, Inc.



4.0mm Fasteners

Part Numbers

340-6014 to 340-6060

Lengths

14-60mm, 2 mm increments

4.0mm Fastener Instruments

Part Numbers

110103
120005
120007

Description

2.5mm Hex Driver w/AO Adapter
2.8mm Drill 165mm, Calibrated
4.0 x 110mm Standard Tap

3.5mm Fasteners

Part Numbers

335-1010 to 335-1050
335-1055
335-1060

Lengths

10-50mm, 2 mm increments
55mm
60mm

3.5mm Fastener Instruments

Part Numbers

110013
110014
110028
110444
110103

Description

2.5mm Drill Bit 110mm
2.5 x 145mm Drill Bit
3.5mm Drill Bit
3.5mm Tap
2.5mm Hex Driver w/AO Adapter

Proximal Humerus Plates

Part Numbers

316-03-090
316-05-115
316-07-140

Number of Holes

3
5
7

Length (mm)

90mm
115mm
140mm

Other Instruments

Part Numbers

110205
110302
110310
110500
110522
110702
120001
120002
120003
120004
120006
120008
120009
120010
120011
130010

Description

Depth Gauge, 60mm
Silicone Handle, AO QC, 25mm x 135mm
QC T-Handle
Countersink
3.5mm/2.5mm Universal Drill Sleeve
2.5mm Hex Driver with Silicone Handle
K-Wire Sleeve
Drill Sleeve
Screw Insertion Sleeve
Wire Depth Gauge
Threaded Drill Guide
Metaphyseal Drill Guide
1.5Nm Torque Limiter
3.5mm Hex Driver
Metaphyseal Drill Guide, Short
1.60, K-wire, Trocar, Smooth, 150mm

References:

DeBaun, M., Swinford, S., Chen, M., Thio, T., Behn, A., Lucas, J., Bishop, J. and Gardner, M., 2020. Biomechanical comparison of bone-screw-fasteners versus traditional locked screws in plating female geriatric bone. Injury, 51(2), pp.193-198.

Ha Vo, Lawrence Webb, Bich Nguyen, Trung Le, "SMV Orthopedic Company Bone-Screw-Fastener Resistance Against Multi-Directional Forces and Bending Moments", Orthopaedic & Tissue Mechanics Laboratory, Mercer University, Date of the Study.



OsteoCentric®
TECHNOLOGIES



UnifiMI™
An OsteoCentric
Technology

75 West 300 N, Suite 150
Logan UT, 84321
Phone: 1-800-969-0639
info@osteocentric.com
osteocentric.com

OsteoCentric Trauma, LLC; OsteoCentric SI Fusion, LLC; OsteoCentric Extremities, LLC; OsteoCentric Spine, LLC; OsteoCentric Sports LLC; OsteoCentric Recon, LLC; OsteoCentric Dental, LLC; OsteoCentric Oncology, LLC; and OsteoCentric Vet, LLC; are a family of the companies under the OsteoCentric brand and are under common ownership or common corporate control within OsteoCentric Technologies, Inc.

TR-008-LBL-0002 Rev 05