

OsteoCentric Technologies

Headless Compression System



Discover Unparalleled Stability via Integrated Fixation

- Instant primary stability created by innovative thread geometry and cutting flute technology.
- Consistent major and minor diameter to provide seamless intra-operative removal and adjustments.
- Variable pitch interlocking threads designed to optimize compression and resist off-axis loads.



2.8mm Micro Screw Fasteners

Part Numbers

EX01-5028-12-30

Length

12-30mm, 2 mm increments

3.9mm Mini Screw Fasteners

Part Numbers

EX01-5039-16-34

EX01-5039-40-50

Length

16-34mm, 2 mm increments

40-50mm, 5 mm increments

Instruments

K-Wires

Description

Size

EX01-0026

Smooth

0.9mm

EX01-0028

Smooth

1.1mm

EX01-0029

20mm threaded

1.1mm

Drill Bits & Driver Shafts

Description

Size

EX01-0030

Cannulated Drill for 2.8mm

2.1mm

EX01-0031

Cannulated Drill for 3.9mm

2.9mm

EX01-0032

Cannulated Driver for 2.8mm

1.5mm Hex

EX01-0033

Cannulated Driver for 3.9mm

2.0mm Hex

Parallel Wire Guides

Description

Size

EX01-0015

Center to Center Distance

4mm

EX01-0016


Center to Center Distance

5mm

EX01-0021

Center to Center Distance

6-7-8mm



75 West 300 N, Suite 150
Logan UT, 84321
Phone: 1-800-969-0639
info@osteocentric.com
osteocentric.com