


OsteoCentric®

Pedicle Fastener System



Featuring:



UnifiMI™
An OsteoCentric
Technology

Mechanical Advantages of UnifiMI® Technology

UnifiMI Technology in a pedicle fastener application allows for an increase in stiffness in all planes including axial, lateral medial bending, flexion, and extension compared to the standard legacy thread option.¹



OsteoCentric
SPINE

OsteoCentric® Pedicle Fastener System

UnifiMI® Interface Technology Enabled Implants




Design Rationale & System Description

- Removable constraint collar provides extra control
- Visual guidance for rod delivery and construct locking
- Reliable, over-the-skin rod length measurement
- Common implant interfaces to reduce instrumentation
- 30mm integrated rod reduction capability
- Robust hybrid square thread minimizes cross threading during rod reduction

Enhanced by UnifiMI®

UnifiMI creates a 58% increase in implant off-axis stability².



Unique thread geometry instantly and circumferentially interlocks with bone by entrapping and containing bone between the thread form.

This mechanical interlocking creates a structural and functional connection between an implant and bone which performs similar to Biological Integration (Osseointegration).

OsteoCentric Fasteners are designed to create a new standard for primary implant stability and expedite robust long-term stability for implants utilizing Osseointegration (OI).



This document is intended exclusively for physicians and is not intended for laypersons. Information on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendations. Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and/or advice.

All content herein is protected by copyright, trademarks and other intellectual property rights owned by or licensed to OsteoCentric Technologies, Inc. or one of its affiliates and must not be redistributed, duplicated, or disclosed, in whole or in part, without the express written consent of OsteoCentric Technologies, Inc.

Pedicle Fastener Implant List

Part#	Description	Set Quantity
101-10001	Set Screw	14
101-10035	5.5 x 35mm Curved Rod	2
101-10040	5.5 x 40mm Curved Rod	2
101-10045	5.5 x 45mm Curved Rod	2
101-10050	5.5 x 50mm Curved Rod	2
101-10055	5.5 x 55mm Curved Rod	2
101-10060	5.5 x 60mm Curved Rod	2
101-10065	5.5 x 65mm Curved Rod	2
101-10070	5.5 x 70mm Curved Rod	2
101-10075	5.5 x 75mm Curved Rod	2
101-10080	5.5 x 80mm Curved Rod	2
101-10090	5.5 x 90mm Curved Rod	2
101-10100	5.5 x 100mm Curved Rod	2
101-10110	5.5 x 110mm Curved Rod	2
101-10120	5.5 x 120mm Curved Rod	2
101-55035	5.5 x 35mm Cann Extension Fastener	2
101-55040	5.5 x 40mm Cann Extension Fastener	4
101-55045	5.5 x 45mm Cann Extension Fastener	4
101-55050	5.5 x 50mm Cann Extension Fastener	2
101-65035	6.5 x 35mm Cann Extension Fastener	4
101-65040	6.5 x 40mm Cann Extension Fastener	8
101-65045	6.5 x 45mm Cann Extension Fastener	8
101-65050	6.5 x 50mm Cann Extension Fastener	6
101-65055	6.5 x 55mm Cann Extension Fastener	4
101-75040	7.5 x 40mm Cann Extension Fastener	4
101-75045	7.5 x 45mm Cann Extension Fastener	4
101-75050	7.5 x 50mm Cann Extension Fastener	2
101-75055	7.5 x 55mm Cann Extension Fastener	2

* Approved for additional lengths. Please inquire regarding availability *



5.5mm Fastener



6.5mm Fastener



7.5mm Fastener

Pedicle Fastener Instrument List

Part#	Description	Set Quantity
101-0002	10.1 Nm Torque Limiter Handle	1
101-0004	Large Head Positioner	1
101-0006	T25 Hudson Tapered Insert	2
101-0007	Cannulated Probe Handle	1
101-0008	Cannulated Probe Shaft	1
101-0012	Cann Tap for 5.5mm Fastener	1
101-0013	Cann Tap for 6.5mm Fastener	1
101-0014	Cann Tap for 7.5mm Fastener	1
101-0016	Bullet Rod Introducer	2
101-0017	Bullet Rod Introducer Shaft	2
101-0019	Extension Counter Torque	1
101-0021	Rod Length Template	1
101-0022	T25 Driver, Square	1
101-0023	Tab Removal Tool	1
101-0024	T25 Driver, Hudson	1
101-0025	T25 Driver, 1/4 Square	2
101-0026	Compressor Top Handle	1
101-0027	Compressor Bottom Handle	1
101-0028	Distractor Top Handle	1
101-0029	Distractor Bottom Handle	1
101-0030	Restraint Tube	6
101-0031	MIS T25 Super Locking Driver	2
101-0032	Straight Ratcheting Handle	2
101-0033	Open T25 Super Locking Driver	2
101-0034	Small Pedicle Sounder	1
101-0035	Large Pedicle Sounder	1
101-0036	Curved Lenke Probe	1
101-0037	Tab Breaker Block	*
101-0038	Rescue Sleeve	*
101-0112	Solid Tap for 5.5mm Fastener	*
101-0113	Solid Tap for 6.5mm Fastener	*
101-0114	Solid Tap for 7.5mm Fastener	*
101-90000	Sterilization Case Lid	2
101-90001	Rod Module Lid	1
101-90003	Rod Module	1
101-90004	MIS Implant Case	1
101-90005	MIS Instrument Case	1
101-90007	MIS Implant Case Upper Tray	1
101-90008	MIS Implant Case Lower Tray	1
101-90009	MIS Instrument Case Upper Tray	1
101-90010	MIS Instrument Case Lower Tray	1
104-0006	Small Hudson Power Adapter	*
201-0009	Small Dilator Tube	2
201-0010	Medium Dilator Tube	2
201-0011	Large Dilator Tube	2
909-0480N	480mm Nitinol Guidewire	16

*These instruments are only available upon request and are not present in the system case and tray.

Reference:

1. Nguyen B, Le T, Vo H, Speece B, Webb L. A Novel Method for Evaluating Pedicle Screws Utilizing UnifiMI Thread Technology. Mercer Biomedical Engineering Society. Biomechanical Testing Report.
2. OsteoCentric testing data on file and numbers reflect specific tests in various medium and are not meant to reflect generic claims for improved performance across all clinical applications. Approved conditionally and non-product specific.



75 West 300 N, Suite 150
Logan UT, 84321
Phone: 1-800-969-0639
info@osteocentric.com
osteocentric.com

OsteoCentric Trauma, LLC; OsteoCentric SI Fusion, LLC; OsteoCentric Extremities, LLC; OsteoCentric Spine, LLC; OsteoCentric Sports LLC; OsteoCentric Recon, LLC; OsteoCentric Dental, LLC; OsteoCentric Oncology, LLC; and OsteoCentric Vet, LLC; are a family of the companies under the OsteoCentric brand and are under common ownership or common corporate control within OsteoCentric Technologies, Inc.