


OsteoCentric® Trauma

Advanced Pelvic Fixation System



OsteoCentric®
TECHNOLOGIES

Advanced Pelvic Fixation Technique System Description (APFT)

The Advanced Pelvic Fixation Technique System was designed to provide additional intraoperative optionality when reducing and stabilizing complex fractures. These additional instruments complement the current OsteoCentric 7.0mm/8.0mm Large Cannulated Screw Fastener System.

Benefits

- Provides additional instrumentation to aid in fracture reduction and fixation
- Features a threaded and non-threaded beveled tip guide wire for optimal maneuverability
- Features low profile drill sleeve for a minimally invasive approach




APFT Instruments

Instruments are made from various materials approved per ASTM F899.

Part Number	Description
2004-0101	2.7x250mm Transsacral Beveled Tip Wire
2004-0102	2.7x350mm Transsacral Beveled Tip Wire
2004-0103	2.7x450mm Transsacral Beveled Tip Wire
2004-0141	2.7x250mm Transsacral Beveled Tip Wire THD
2004-0142	2.7x350mm Transsacral Beveled Tip Wire THD
2004-0143	2.7x450mm Transsacral Beveled Tip Wire THD
2004-0111	4.5x270mm Drill Bit, Side Cutting
2004-0121	4.5mm Drill Sleeve, 100mm
2004-0122	4.5mm Drill Sleeve, 160mm
2004-0131	Mallet
2004-0023	2.7mm Thd Toracr Tip Guidewire
202-80100	2.0mm Guide Wire 254mm
320-80032*	1.6mm Guide Wire Trocar 230mm
2004-2000-17	Advanced Pelvic Fixation Technique Caddy

Available In Threaded & Non-threaded Options



**These instruments are only available upon request and are not present in the system case and tray.*

This document is intended exclusively for physicians and is not intended for laypersons. Information on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendations. Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and/or advice.

All content herein is protected by copyright, trademarks and other intellectual property rights owned by or licensed to OsteoCentric Technologies, Inc. or one of its affiliates and must not be redistributed, duplicated, or disclosed, in whole or in part, without the express written consent of OsteoCentric Technologies, Inc.



75 West 300 N, Suite 150
Logan UT, 84321
Phone: 1-800-969-0639
info@osteocentric.com
osteocentric.com

OsteoCentric Trauma, LLC; OsteoCentric SI Fusion, LLC; OsteoCentric Extremities, LLC; OsteoCentric Spine, LLC; OsteoCentric Sports LLC; OsteoCentric Recon, LLC; OsteoCentric Dental, LLC; OsteoCentric Oncology, LLC; and OsteoCentric Vet, LLC; are a family of the companies under the OsteoCentric brand and are under common ownership or common corporate control within OsteoCentric Technologies, Inc.