

OsteoCentric

UnifiMI 3.5mm Long Screw Fastener Set


Superior Stability Performance in Multi-axial Load Environments

First ever long fastener featuring UnifiMI technology to mechanically integrate with the bone to empower surgeons with instant primary stability.

UnifiMI technology resists multiplanar and multidirectional loads to create Integrated Fixation.

Maintains purchase strength during and after cyclic loading.

Reduced toggle and implant displacement.




OsteoCentric
TECHNOLOGIES

Reference: Ha Vo, Lawrence Webb, Bich Nguyen, Trung Le, "SMV Orthopedic Company Bone-Screw-Fastener Resistance Against Multi-Directional Forces and Bending Moments," Orthopaedic & Tissue Mechanics Laboratory, Mercer University, May 2019.

OsteoCentric 3.5mm Long Screw Fastener Set featuring UnifiMI Technology

3.5mm Long Screw Fasteners

Part Number	Length
335-1052-88	52-88mm, 2 mm increments
335-1090-1130	90-130mm, 5 mm increments



Instruments

110028	3.5mm Drill Bit, 110mm
110040	2.5mm Drill Bit, 200mm, Calibrated
110104	Long 2.5mm Hex Driver AO
110206	Depth Gauge, 150mm
110454	3.5mm Cortical Screw Tap, 200mm
110701	2.5mm Hex Driver Sleeve
110703	Long 2.5mm Hex Driver, Fixed Silicone Handle 300mm
110705	SCREW HOLDING FORCEPS



75 West 300 N, Suite 150
Logan UT, 84321
Phone: 1-800-969-0639
info@osteocentric.com
osteocentric.com